





Technische Information, Technical Information, Informations techniques,  
Información técnica, Informazioni tecniche, Teknisk information, Техническа  
информация, Technische informatie, Informacje techniczne, Τεχνικές  
Πληροφορίες, Technické informace, Műszaki információk, Teknik Bilgiler,  
Техническая информация



**MD**

**4W9... / 4W8... / 4W70.. / 4L9...**

DE - LITE Liner® Silikon ohne Stoffbezug	4
GB - LITE Liner® Gel with fabric cover	11
FR - LITE Liner® Silicona con funda de tela	18
ES - LITE Liner® Gel con funda de tela	25
IT - LITE Liner® Gel con rivestimento in tessuto	32
SV - LITE Liner® Gel med tygöverdrag	39
BG - LITE Liner® Gel с покритие от плат	46
NL - LITE Liner® Gel met stoffen bekleding	54
PL - LITE Liner® Gel z osłoną z tkaniny	61
EL - LITE Liner® Gel με υφασμάτινο κάλυμμα	68
CS - LITE Liner® Gel s látkovým potahem	76
HU - LITE Liner® gél szövet borítással	83
TR - Kumaş kaplamalı LITE Liner® Jel	89
RU - LITE Liner® Gel с тканевым покрытием	96

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MD

### Informationen

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## GB LITE Liner® Gel with fabric cover

#### Information

- Read this document carefully before using the product. Follow the safety instructions to avoid injury and product damage.
- Instruct the user in the proper and safe use of the product.
- Reporting a serious incident: Any serious incident related to the product must be reported to the manufacturer and the competent authorities.
- Keep this document in a safe place.

#### 1. Description

The product is a liner made of a theroplastic TPE gel that covers the residual limb and provides a contact surface between the residual limb and the socket of the prosthesis. It is available as a locking liner with distal connection option or as a cushion liner without distal connection option.

A locking liner has either a fixed or a removable distal connector.

A cushion liner does not have a distal attachment and must be used with a sleeve to hold the vacuum.

#### 2. Indicated Use

The Gel liners are to be used exclusively for the prosthetic fitting of the lower extremities. The product is the interface of the residual limb with the prosthetic socket.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

#### It fulfills the following functions:

- The prosthesis is attached to the limb stump using pin locking or vacuum system technology
- residual limb compression (prevention of edema, reduction of residual limb volume fluctuations)
- Protection and care of the skin
- Shock absorption
- Cushioning

#### 2.1 Indication

The Gel liners are only to be used for prosthetic restoration after amputation and/or congenital deficiency of the lower extremities.

- Poor soft tissue coverage
- Prominent bone protrusions

soft tissue substitution by the elastomer

The gel flows around such pressure zones and increases the comfort of the residual limb bedding

- bums (skin grafts)

The liner absorbs a large part of the shear forces when walking and thus largely protects these sensitive skin zones from mechanical overload.

- Skin adhesions on the periosteum

shear forces are kept away from the corresponding skin zones to a certain extent.

## 2.2 Contraindication

- Too many soft tissues/soft tissue overhangs  
The patient „swims“ in his stump und shows an unsafe walking.
- Unstable knee joint  
All forms of short prostheses clearly produce instabilities of the knee joint and should be applied only with medical and physiotherapeutic support. Such patients will generally complain about a too high axis rise and a low rotational stability between socket and stump. Furthermore, the durability of the liner will be significantly impaired. Such patients should be fitted with a matrix liner.
- Users with a high activity level  
If a stump end bearing is not possible, the LITE Liner® Gel must not be used as the stump end must normally bear 20-30 % of the body weight. Please follow the instructions for cleaning
- No residual limb end loading possible
- Bad hygiene
- Very sensitive neuromas
- Allergy / intolerance to the material
- Contact with irritated / sick skin
- Contact with open wounds

If you have severe diabetes and poor blood circulation, consult a doctor before use

## 3. Technical data

Name	Art. No.	Activity Level	TT / TF
LITE Liner® GEL SUPREME Locking	4W9000, 4W9100	1-3	TT
LITE Liner® GEL SUPREME Matrix Locking	4W9010, 4W9110	1-3	TT
LITE Liner® GEL SUPREME Cushion	4W9500, 4W9600	1-3	TT
LITE Liner® GEL SUPREME Matrix Cushion	4W9510, 4W9610	1-3	TT
LITE Liner® GEL SUPREME TF Locking	4W90TF, 4W91TF	1-3	TF
LITE Liner® GEL SUPREME TF Cushion	4W95TF, 4W96TF	1-3	TF
LITE Liner® GEL SUPREME TT Locking	4W90TT	1-3	TT
LITE Liner® GEL SUPREME TT Cushion	4W95TT	1-3	TT
LITE Liner® GEL PRIME Locking	4W8000, 4W8100	1-3	TT
LITE Liner® GEL PRIME Matrix Locking	4W8010, 4W8110	1-3	TT
LITE Liner® GEL PRIME Cushion	4W8500, 4W8600	1-3	TT
LITE Liner® GEL PRIME Matrix Cushion	4W8510, 4W8610	1-3	TT
LITE Liner® GEL EXCELSIOR Locking	4W7000	2-4	TT
LITE Liner® GEL EXCELSIOR Cushion	4W7500	2-4	TT
LITE Liner® GEL EXCELSIOR TF Locking	4W70TF	2-4	TF
LITE Liner® GEL EXCELSIOR TF Cushion	4W7500	2-4	TF
LITE Liner® GEL GELACTIVE Locking	4L9010	2-4	TT
LITE Liner® GEL GELACTIVE Cushion	4L9510	2-4	TT

### 3.1 Liners with distal connector

These liners can be used for nearly every stump shape. We recommend a total surface weight bearing socket with distal end contact in combination with a locking system.

Patients with a very high level of activity will generally complain about a too high axis rise and a low rotational stability between socket and stump. Such patients should be fitted with a matrix liner.

### 3.2 Liners without distal connector

These liners can be used for nearly every stump shape. We recommend a total surface weight bearing socket with distal end contact in combination with an air outlet valve.

### 3.3 Service life

The product is a wearing part that is subject to normal wear and tear. After a period of use of 6 months, we recommend replacing the liner.

### 3.4 Restriction of Use

Due to the adjustment, which must be carried out by trained specialists, and for reasons of hygiene, we would like to point out that this product is only to be used by one person.

### 3.5 Misuse

The product is intended exclusively for prosthetic fitting in accordance with its intended use and is not suitable for other sports activities.

### 3.6 Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- Only combine the product with prosthetic components that are approved for that purpose.
- Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

### 3.7 Use under unallowable environmental conditions

Risk of injury due to damage to the product

- Do not expose the product to unallowable environmental conditions.
- If the product has been exposed to unallowable environmental conditions, check it for damage.
- If damage is apparent or in case of doubt, do not continue using the product.
- Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

### 3.8 Environmental conditions

Storage and transportation

Do not expose the product to temperatures above +55 °C and below -20 °C.

Items must be stored at room temperature in a dry place out of direct sunlight if left unused for a long time.

## 4. General safety instructions

Keep the product safe from children - There is a risk of suffocation if the product is pulled over the mouth and nose!

Take care not to wear the liner on injured skin!

The Gel liner should not be worn directly on open wounds or broken skin, or used by people prone to skin inflammation.

Upon receipt and at regular intervals, please check the edge of your gel liner for any tears or irregularities. Such defects in the edge area can lead to cracks in the material and damage to your TPE gel socket. Improper use or use of a defective product can lead to potential adhesion failure. If such defects occur, the user must stop using the product and must consult an orthopedic technician.

To avoid the risk of personal injury and product damage, do not use the product beyond the tested service life.

Incorrect selection, adjustment, use and/or inadequate control can lead to damage to health.

The prosthesis and its functional components must be checked regularly by an orthopedic technician for function, possible wear and tear and possible damage.

If components for different activity levels are combined, the activity level of the component with the lowest activity applies to the entire prosthesis.

## 5. Size selection

### 5.1 How to choose the correct liner size

- Measure the circumference of the residual limb at 4 cm from the distal end (see figure 1). With a transtibial liner, the knee should be slightly bent to allow the soft tissue to hang down.
- Choose the liner size corresponding to this measurement or the next smaller size, respectively.
- Do not choose a liner whose size exceeds the measurement taken.
- Never reduce by more than one liner size in order to avoid unnecessarily high compression.

A liner with a connector often has a matrix with a length of 10 cm integrated. This matrix should end below the head of the fibula and must not exceed it, as this could result in problems during knee flexion.

The selection of the correct liner size is of crucial importance for comfort and safety.

- If the liner is too tight, this could result in pistoning, numbness and blisters.
- If the liner is too loose, this may cause increased perspiration and rotation of the stump inside the liner which could lead to blisters and skin rashes.
- Please contact your prosthetist if one of the symptoms or other problems should occur.

## 6. Adjustment of the liner

Threaded liners usually have a matrix. The matrix prevents excessive elongation of the liner, with little restriction of radial extensibility. The length of the matrix depends on the application area of the liner. Liners with a matrix in the distal area are particularly suitable for TT prostheses. Liners with a full-length matrix are commonly used for TF prostheses.

The product must be thoroughly inspected for defects before fitting.

### 6.1 Trimming of a liner

Do not follow socket edge during trimming (see figure 2), but always cut all around at an equal height slightly above/beyond the socket edge. Do not use zigzack scissors in order to avoid a possible tearing of the heavily pre-stressed liner at the serrated areas during trimming. Do not trim the liner below the edge of the matrix or below the proximal stem cut lines. Excessive trimming may reduce adhesion between the product and the limb, thereby affecting fixation.

### 6.2 Customization of the liner

Thermoformability of the copolymer liner (see figure 3)  
Roll the liner onto a clean and preheated plaster positive. Then put it in an oven for about 30 minutes at 65-70°C. Afterwards, allow the liner to cool down on the plaster positive and remove it from the positive once it has completely cooled down.

Only the positive model should be in contact with the oven. A direct contact between the LITE Liner® Gel and any part of the oven may result in a damage of the LITE Liner® Gel. Lengthen your positive model about 10-15 cm at the proximal end. By means of this you create a support area and avoid damages on the gel.

### 6.3 Adjustment of the socket

Proper functioning of the product depends to a large extent on the fit of the upper. The following points should be checked in order to achieve the highest level of adhesion, control and comfort.

- The patient can easily put on the socket without the use of lubricants.
- Complete contact between the liner and the inner wall of the socket.
- No development of air pockets when the patient puts on the prosthesis.
- Allows full movement without compromising distal adhesion.

If user feedback indicates that one or more of the above conditions are not met, the socket must be redesigned or remanufactured accordingly

When using a cushion liner, the socket must be airtight to ensure vacuum fixation when a sleeve is applied.

### 6.4 Optional: Installing the pin

The pin for connecting the liner to the shuttle lock is screwed into the thread at the distal end of the liner. Apply medium-strength Loctite to the thread beforehand. Screw the pin into the thread and tighten carefully with 3-4 Nm.

At a fitting with a locking system, please make sure that the patient has a socket with distal end contact. Always position the locking system in line with the residual limb center.

## 7. Description

### 7.1 Donning of the liner

Cover open wounds or broken skin with a plaster or other appropriate wound dressing to prevent direct contact between the wound and the liner.

1. Turn the gel liner inside out, expose the distal end of the liner as much as possible and place it against the stump by avoiding any air pockets (see figure 4).
2. Take care that the screw thread is centered and in line with the stump (see figure 5).
3. Then roll the liner on the stump carefully (see figure 6). To avoid unpleasant wrinkles in the area of the hollow of the knee, carefully pull the socket up behind the knee and put it on the skin again in a stretched state. Take special care that the socket will not be damaged and that the pin receptacle is positioned as central as possible.
4. Please avoid to pull at the upper end of the liner (see figure 7) as this results in an increased tension on the liner brim and the subjacent skin. This manifests itself in water blisters, weals or similar symptoms. If wrinkles should have occurred during donning of the liner, please repeat procedure.

- When donning the LITE Liner® Gel at the first time, the patient may feel it as slightly tight. To achieve a more comfortable fit, please proceed as follows:  
Roll the LITE Liner® Gel, after having folded over the upper edge, nearly completely down the stump so that the liner looks like a thick loop. Stop rolling down just before the distal end and subsequently roll the liner up the stump again. If the LITE Liner® Gel is, then, still too tight, an additional customization (thermoplastic deformation) may be necessary.
- Avoid wrinkles and trapped air between the stump and the LITE Liner® Gel.
- In case there is still space between the LITE Liner® Gel and the socket, you should either use a thicker LITE Liner® Gel or a stump sock or an additional soft sock.
- Then let the patient step into the socket. If using a liner with distal connector, please observe the corresponding instructions for the locking mechanism used by your patient.

In the case of TPE gel liners, it is advisable to treat the skin in areas of high friction with skin protection lotion 4W296 or selectively with Protect-Gel 4W293 before putting on the liner.

A large-area treatment with Protect-Gel 4W293 may cause a losing of the liner. The application of Protect-Gel 4W293 is intended only for small areas, i.e. as big as a coin. We, therefore, implicitly recommend to apply the corresponding care products at least 10 minutes before donning of the liner.

### 7.2 Perspiration

During the first weeks of wearing a LITE Liner® Gel it may happen that the patients complain about an increased perspiration. Unscented antiperspirants 4W297 or alum which are applied onto the stump before donning the LITE Liner® Gel, are moderate and effective means to avoid excessive perspiration. Please observe the corresponding instructions on the product leaflet and never use these products on open wounds on the stump.

### 7.3 Donning the socket

#### Cushion liner:

Use this in conjunction with a knee cap. Refer to the knee cap user manual for instructions on correct use.

Using your body weight, slide into the socket until the residual limb is firmly seated before using the knee cap to establish fixation in the Vacuum System technique.

#### Locking liner:

Using body weight, slide into the shaft until the distal connector firmly engages the detent and a secure connection is made.

Make sure the connection between the distal connector and the detent is secure and free of play.

## 8. Cleaning and care

### 8.1 Care of the liner

Clean the product daily after use and before first use.

Liners that have already been worn usually feel softer than new liners.

A proper hygiene is extremely important when using the LITE Liner® Gel. The LITE Liner® Gel must be thoroughly cleaned in the evening after wearing and should not be worn for the next 24 hours so that the perspiration which is possibly in the material, can ventilate. For this reason, it is advisable, to fit every patient with two LITE Liner® Gel.

### 8.2 Daily cleaning

- Roll the LITE Liner® Gel off the stump so that the gel side is then facing out.
- Clean the LITE Liner® Gel immediately after removing – with the gel side still facing out – with lukewarm tap water and a ph-neutral soap. Apply the soap by hand or with a clean and soft cloth or sponge (see figure 8&9).

Do not treat the LITE Liner® Gel with a brush. This could roughen the surface of the liner which may lead to skin irritations.

- Please see that, afterwards, all soap residues are rinsed from the LITE Liner® Gel.
- Towel the gel side of the LITE Liner® Gel with a clean and lint free cloth (see figure 10).
- Turn the LITE Liner® Gel inside out so that the fabric is pointing outwards and the gel side is pointing inwards.
- Put the LITE Liner® Gel - with the fabric side pointing outwards - over the drying stand.
- Then clean the stump.

### 8.3 Weekly cleaning (for disinfection)

1. Apply a small amount of ethyl- or isopropyl alcohol on a soft and clean cloth.
2. Now wipe the gel side of the LITE Liner® Gel with the alcohol-based cloth for about 2 minutes.
3. Put the LITE Liner® Gel – with the fabric side pointing outwards- over the drying stand.

Do not soak the LITE Liner® Gel in ethyl- or isopropyl alcohol as such an extended contact with ethyl- or isopropyl alcohol can lead to a stiffening of the liner.

Using the wrong detergent may damage the product.

Do not use bleach, harsh chemicals or objects such as salt water, acids, dust or sand as this may damage the fabric on the outside of the liner or the gel surface.

After each brief exposure to salt or chlorinated water, the inside and outside of the LITE Liner® GEL must be thoroughly rinsed with fresh water. The LITE Liner® GEL must then be dried.

Dry the liner textile side out to prevent damage.

When handling your gel socket, be careful not to damage it with your fingernails or sharp objects such as the locking pin.

Please also avoid contact of the liner with glass or carbon fibers or other foreign particles, as these materials cannot be removed by normal cleaning. Any residue left on the gel can cause skin irritation. If the liner has come into contact with such materials, it should be returned to the prosthetist.

### 8.4 Storage of the LITE Liner® Gel

If the LITE Liner® Gel is not worn, it should be stored at a cool and dry place. If the LITE Liner® Gel is not used for a longer period of time, so it is recommended to cover it with a clean plastic bag.

### 8.5 Care of the stump

At the beginning of your fitting with a gel liner, you may experience increased perspiration. However, this will decrease significantly after one to two weeks, since, in contrast to more conventional fittings, there is hardly any friction between the skin and the socket. With increased physical exertion, sweating can still occur. In the case of sensitive skin, skin irritation can occasionally occur. In this case, consult the medical supply store that carried out your treatment or your doctor. Such skin irritation is caused, for example, by the slight pull on the hair, which can irritate the hair roots.

The stump must be thoroughly washed daily with neutral soap. We recommend using a mild, pH-neutral, fragrance-free and dye-free soap. For dry skin, apply a pH-balanced, fragrance-free and dye-free lotion to condition and soften the skin.

## 9. Additional products

The following accessories can be used in combination with cushion and locking liners.

### 9.1 Use of Distal Cups 4W2005

The Distal Cups 4W2005 offer increased comfort for conical stumps with small soft tissue covering. They are used in combination with gel liners in order to compensate conical stump shapes and provide an additional padding.

The Gel Distal Cups 4W2005 adapt to the stump and, consequently, decrease pressure over sensible bony areas during walking.

### 9.2 Use in combination with the Volume Management Pad 4W17

If there are cavities in the distal area of the liner, e.g. due to scar indentations in the residual limb or short-term weight fluctuations in this area, these can be compensated for using one or more volume management pads 4W17. To do this, the gel pad corresponding to the size of the liner is selected and positioned in the areas of volume loss before the liner is put on.

The liner is worn over the volume pad.

1. position the volume pad on the residual limb with the thin, beveled side facing the skin.
2. if necessary, volume pads of graduated sizes can be used on top of each other to compensate for larger volume fluctuations.

Please ask your prosthetist to demonstrate this procedure exactly in order to avoid an improper or careless application which may result in skin irritations or even in hematomas.

### 9.3 Use of stump socks

Any fluctuations in volume of the residual limb that may occur should be compensated for if there is a noticeable loosening of the seat in the prosthesis socket.

Stump socks with a hole for the pin are available for liner variants with a PIN connection.

For cushion liner variants, we recommend stump socks without a distal hole.

To protect the gel liner from possible damage caused by friction, a thin residual limb stocking with a distal opening for the pin should also be pulled over the gel liner. This reduces the friction between the socket and the liner.

#### 10. Maintenance/Service

The Gel liner is subject to natural wear and tear during continuous use and must therefore be checked at regular intervals and replaced if necessary.

##### **CAUTION: Failure to follow maintenance instructions**

Risk of injury due to a change in function or loss of function as well as damage to the product  
Observe the following maintenance instructions:

- Perform a visual and functional check of the prosthesis components after the first 30 days of use.
- Check the complete prosthesis for wear during the normal consultation.
- Carry out annual safety checks.

#### 11. Guarantee

The following guarantee declaration applies in addition to any existing legal warranty claims and is independent of these.

##### 11.1 Conditions of entitlement

The manufacturer provides a warranty on the product from the date of purchase.

A guarantee can only be given if the product is used under the specified conditions and for the intended purposes. The warranty covers defects that are demonstrably due to material, manufacturing or construction errors and that are claimed against the manufacturer within the warranty period. Claims under the manufacturer's warranty can only be made if the product concerned has been installed or adapted by or under the supervision of trained specialist personnel. Maintenance intervals, if any, must have been demonstrably carried out.

##### 11.2 Scope and exclusion of warranty

The product is guaranteed to be free of material, manufacturing or construction defects. The manufacturer's warranty does not cover defects that are attributable :

- to incorrect installation
- on operating errors
- to improper / improper use
- on normal wear and tear
- on improperly machined or sharp-edged edges of the socket
- external influences such as force majeure or third party intervention

##### 11.3 Duration of warranty

The duration of the manufacturer's warranty varies depending on the product. In any case, the duration of the manufacturer's warranty starts from the date of delivery to the customer. If used properly, the fabric coated LITE LINER Gel is guaranteed for 6 months.

##### 11.4 Liability

The manufacturer is only liable if the product is used in accordance with the descriptions and instructions in this document. The manufacturer is not liable for damage caused by non-compliance with this document, in particular by improper use or unauthorized modification of the product. Each manufacturer can only be held liable for the failure of its own fitting parts. A further liability of the manufacturer is only possible if his own components can be proven to have been the cause of the damage or functional failure of components from other manufacturers.

##### 11.5 Warranty service

If all conditions for a warranty claim are met, Wagner Polymertechnik GmbH will, at its discretion, repair the damaged product, send a replacement free of charge or issue a credit note for the cost of the product.

#### 12. Declaration of Conformity

As the manufacturer, Wagner Polymertechnik GmbH declares under sole responsibility that the LITE Liner® Gel with fabric cover comply with the requirements of Regulation (EU) 2017/745 for medical devices (MDR).

#### 13. Disposal

The product is to be recycled or disposed of without risk to others or the environment. Please observe the applicable national regulations.

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