Purpose of the instructions

This Manual is relative to a 4-points knee orthosis, in compliance with Regulation 2017/745 concerning class 1 medical devices, and intends to provide guidance for the application, regulation and the proper use of the device. The declaration of conformity was therefore issued by Tecnoway under its responsibility, in accordance with Annex IV of this Regulation.

Identification of the manufacturer and the product

Tecnoway produces medical devices, in particular orthopaedic braces. The Firm is constantly searching for new materials innovative technologies and ergonomic solutions. Its system of quality management is certificated in accordance with ISO 9001 and ISO 13485 regulation. The product specified in this manual is a knee orthosis used for the post-operative treatment of the ligaments reconstruction of the knee, as limitation of the knee-joint in the post-traumatic rehabilitation of the knee.

General description and mode of action

The device described in this manual is composed primarily of a frame in aluminum alloy and non-toxic material padded supports making the device light, comfortable and safe. Easily removable, it controls the movements of the knee without restricting the operations of daily hygiene

Caution for use

The device should be used exclusively for orthopaedic treatment of the knee. The duration of daily use and the application period must be indicated in the prescription of an orthopaedic doctor. The pressure exerted by the device should not act upon injury, or caused localized symptoms caused by compression, or compress nervs fibers and blood vessels. The device is not fireproof, do not use near flames.

Features

Anatomical structure frame in light alloy of aluminium anticorodal, coated by MTP padded with a soft-grip fabric inside to reduce perspiration and the migration of the orthosis. Separate bands in order to ensure greater breathability of the leg. Closure Velcro straps by double numbered velour to facilitate the application of the patient. Joint adjustment flexion and extension: flexion 0°/ 10°/ 20°/ 30°/ 45°/ 60°/ 75°/ 90°/ 120°, extension at 0°/ 10°/ 20°/ 30°/ 45°/ 120°. Option to block joint at 0°/ 10°/ 20°/ 30°/ 45°. Condylar protective detachable pads. Its shape provides an excellent anatomic stabilization of the limb and a good tolerability during use.

Indications

Post-operative knee ligaments reconstruction treatment, as limitation of the knee-joint, in the post-traumatic rehabilitation of the knee, such as prevention of sport injuries.

Application (Orthotist)

1. Open the straps to free the VELCRO belts.
2. Position the joints near the knee pivot point so that they are supported by the condyles, making sure that the joints are aligned with each other. (fig A) If the condyloid supports do not adhere to the knee, place the 5 mm thicknesses included (fig. D) between the plastic shell and the condyloid support.
3. Fit the padded straps, tightening them around the leg with the Velcro by following the numerical order. (fig B)
4. Regulate flexion and extension with the desired angulation by using the special metal inserts. (fig C)
5. If necessary adjust manually or with help of tools the metal bars to the morphology of the patient's limb.

Application (for patient)

1. Open the straps to free the velcro belts.
2. Position the joints near the knee pivot point so that they are supported by the condyles, making sure that the joints are aligned with each other. (fig A)

3. Close the Velcro belts and straps. (fig B)

Warning

The device in question must be prescribed by an orthopaedic specialist and applied by an orthotist who will according to the needs and conditions of patient instruct on the proper and safe use of the device. Any modification or regulation different from thre previously set on the device at the time if application and adjustment should be prescribed by an orthopaedic specialist and performed by an orthotist. Don't change the regulation done by the orthopaedic specialist. For hygienic and functional reason we recommend the use for single patient because it could cause risks for hygienic or functional aspect. In particular sensible person, direct contact of skin to device may cause redness and/or irritations, if those problems may occure pls. contact directly your doctor. If the device is no more under normal use it should be replaced and disposed off.

Incident reporting notice

Any serious incident that occurs in connection with this medical device must be reported to the manufacturer and the competent authority.

Storage instructions

When not using the device, even for short periods, store it in its original packaging to maintain a good level of cleanliness. Keep it in a cool, dry place and protect from direct sunlight.

Maintenance

Before use always check the tightening conditions of the screws, of the joints and of the adjustable elements. Before using check the state of wear of the device and its parts. If worn out contact your orthotist to substitute. Hand wash the paddings with warm water and mild soap, rinse accurately. Dry away from heat sources.

Disposal

Make sure that the device is disposed off in a correct manner, do not dispose of in the environment

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This manual has been realized by Tecnoway exclusively for its clients. Its reproduction and/or entirely and/or partial divulgation is forbidden without Tecnoway authorization. This manual comes together with equipment device, and should be kept in an easy access and suitable place. It should be preserved for all the device life cycle. Personal assigned to application, regulation, use and maintenance should carefully read and understand the instructions of this manual. Tecnoway cannot be pursued or held responsible for any damage caused by wrong application different to those in this manual.

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& by **TECNOWAY SRL**
via Montefeltro n. 45 - 61022 Vallefoglia -PU- ITALY
+39.721.496175 info@tecnoway.it www.tecnoway.it
COMPANY WITH SYSTEM CERTIFICATION
ISO 9001 - ISO 13485 - ISO 45001

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